## **Neutron diffraction**




Rob McQueeney














	NEUTRONS		ELECTRONS	
Wavelength range	0.4 - 10 Å	0.1 - 5 Å	0 .04 - 0.2 Å	
Energy range	0.001 - 0.5 eV	3000 - 100000 eV	6000 - 120000 eV	
Cross-section	10 <sup>-25</sup> barns	barns 10 <sup>-25</sup> Z <sup>2</sup> barns ~10 <sup>-22</sup> bar		
Penetration depth	on depth ~ cm ~ μm		~ nm	
Typical flux	10 <sup>11</sup> s <sup>-1</sup> m <sup>-2</sup>	10 <sup>24</sup> s <sup>-1</sup> m <sup>-2</sup>	10 <sup>26</sup> s <sup>-1</sup> m <sup>-2</sup>	
Beam size	mm-cm	μ <b>m-mm</b>	nm-μm	
Typical sample	Any bulk sample	Small crystals, powders, surfaces	Surfaces, thin films, grains, gases	
Techniques	<b>Diffraction</b> Inelastic scattering Reflectivity	Diffraction Photon absorption Photoemission Inelastic scattering	Microscopy Diffraction Emission spectroscopy EELS	
Phenomena	Magnetic/crystal structures collective excitations (phonons, spin waves) electronic excitations (crystal- field, spin-orbit)	Crystal structures, electronic transitions (photoemission, absorption),	microstructure crystal structures electronic transitions	

## **Neutron cross-section**





## Nuclear property


- Random with Z
- Depends on isotope
- Depends on nuclear spin
- Absorption can be problem


		Abundance	Cross-	Absorption
		(%)	section (bn)	(bn)
	Gd		180	49700
	152Gd	0.2	13	735
	154Gd	2.1	13	85
	155Gd	14.8	66	61100
	156Gd	20.6	5	1.5
	157Gd	15.7	1044	259000
	158Gd	24.8	10	2.2
Phy	160Gd	21.8	10.52	0.77
	·	-	· ·	

# Why neutrons??



- Penetration ~ centimeters, bulk probe
- Sensitivity to low-Z
- Isotopic contrast
- Wavelength range ~ interatomic spacing (1-2 Å) → diffraction
  - interacts with nuclei
  - Interacts with magnetic moment of unpaired e-
- Can measure crystal and magnetic structures






# **Producing neutrons**




**Fission**Nuclear reactor



Moderators →Cold-Thermal

# **Spallation**Particle accelerator



Moderators →Cold-Epithermal

Physics 590 6


# **Neutrons by reactor fission**






High flux isotope reactor - ORNL








Physics 590








### High Flux Isotope Reactor at Oak Ridge National Laboratory



ORY

#### The United States' highest flux reactor-based source of neutrons for condensed matter research



07-G00244E/arm

Under consideration





# **Neutrons by pulsed spallation**







## **Spallation Neutron Source (ORNL)**

# **Target-moderator system**





SNS liquid Hg target





Fig. 3. Horizontal cross-section of the flux-trap moderators.



### Spallation Neutron Source at Oak Ridge National Laboratory

#### The world's most intense pulsed, accelerator-based neutron source

#### Backscattering Spectrometer (BASIS) . BL-2

Dynamics of macromolecules, constrained molecular systems, polymers, biology chemistry, materials science

**Spallation Neutrons and Pressure** 

Diffractometer (SNAP) • BL-3

Materials science, geology, earth and

environmental sciences

Chris Tulk + 865.576.7028 + tulkca@ornl.go

Eugene Mamontov · 865.574-5109 · mamontove@ornl.gov

#### Nanoscale-Ordered Materials Diffractometer (NOMAD) • BL-1B (2010)

Liquids, solutions, glasses, polymers, nanocrystalline and partially ordered complex materials

Joerg Neuefeind • 865.241.1635 • neuefeindic@ornl.gov

#### **Wide Angular-Range** Chopper Spectrometer (ARCS) • BL - 18

Atomic-level dynamics in materials science, chemistry, condensed matter sciences

Doug Abernathy · 865.576.5105 · abernathydl@ornl.gov

#### Fine-Resolution Fermi Chopper Spectrometer (SEQUOIA) • BL - 17 (2008)

Dynamics of complex fluids, quantum fluids, magnetism, condensed matter, materials science

Garrett Granroth • 865.576.0900 • granrothge@ornl.gov

#### **Ultra-Small-Angle Neutron Scattering Instrument** (TOF-USANS) • BL-1A (2012\*

Life sciences, polymers, materials science. earth and environmental sciences Michael Agamalian · 865.576.0903 ·

#### Chemical Spectrometer (VISION) • BL-16B (2011)

Vibrational dynamics in molecular systems, chemistry Christoph Wildgruber • 865.574.5378 • wildgrubercu@ornl.gov

BL - 16A

#### **Neutron Spin Echo Spectrometer** (NSE) • BL - 15 (2009)

High-resolution dynamics of slow processes, polymers, biological macromolecules Michael Ohl • 865.574.8426 • ohlme@ornl.gov

#### **Hybrid Spectrometer** (HYSPEC) • BL-14B (2011)

Atomic-level dynamics in single crystals, magnetism, condensed matter sciences

Mark Hagen • 865.241.9782 •

BL-14A

#### Magnetism Reflectometer • BL-4A

Chemistry, magnetism of layered systems and interfaces Valeria Lauter • 865.576.5389 • lauterv@ornl.gov

#### Liquids Reflectometer • BL-4B

Interfaces in complex fluids, polymers, chemistry John Ankner • 865.576.5122 • anknerjf@ornl.gov

#### **Cold Neutron Chopper** Spectrometer (CNCS) • BL - 5

Condensed matter physics, materials science, chemistry, biology, environmental science

Georg Ehlers • 865.576.3511 • ehlersg@ornl.gov

#### **Extended Q-Range Small-Angle Neutron** Scattering Diffractometer (EQ-SANS) • BL-6 (2008)

Life science, polymer and colloidal systems, materials science. earth and environmental sciences

Jinkui Zhao · 864.574.0411 · zhaoj@ornl.gov

**Elastic Diffuse Scattering** Spectrometer (CORELLI) . BL-9 (2013)

> Detailed studies of disorder in crystalline materials Feng Ye • 865.576.0931 • yef1@ornl.gov

**BL-10** 

#### Macromolecular Neutron Diffractometer (MaNDi) • BL-11B (2012)

Atomic-level structures of membrane proteins, drug complexes, DNA

Leighton Coates

#### **Fundamental Neutron** Physics Beam Line • BL-13 (2008)

Fundamental properties of neutrons

Geoffrey Greene • 865.574.8435 • greenegl@ornl.gov

#### **Single-Crystal Diffractometer** [TOPAZ] • BL-12 (2009)

Atomic-level structures in chemistry, biology, earth science, materials science, condensed matter physics

> Christina Hoffmann - 865.576.5127 hoffmanncm@ornl.gov

Powder Diffractometer (POWGEN) . BL-11A (2008)

Atomic-level structures in magnetism, chemistry, materials sciences Jason Hodges • 865.576.7034 • hodgesj@ornl.gov

#### \* Scheduled commissioning date



**Engineering Materials Diffractometer** (VULCAN) • BL-7 (2008)

BL-8A

BL-8B

Mechanical behaviors, materials science, materials processing

Xun-Li Wang • 865.574.9164 • wangxl@ornl.gov










# Places to go






## **Powder diffraction**



For single-crystal:  $Q_{hkl} = k_f - k_i$ 


For powder:  $Q_{hkl} = 2k_i \sin\theta$ 

$$k_i = 2\pi/\lambda$$



Incident beam x-rays or neutrons


Bragg's Law  $\lambda = 2d\sin\Theta$ Powder pattern – scan 2 $\Theta$  or  $\lambda$ 



## **Powder diffraction**




## **Determine the crystal structure**




# **TOF** powder diffraction






POWGEN @ SNS

Time-of-flight  $\tau = L/v = \lambda mL/h = 2mLd\sin\theta/h$ 





Time, wavelength, or d-spacing

# Fitting diffraction data



### Rietveld refinement

- Lattice constants
- space group
- Atom positions
- Site disorder/vacancies
- Thermal vibration amplitude
- Strain broadening

R-factor




Figure 3. Crystal and magnetic structure of the class-I MV (charge-ordered) YBaFe<sub>2</sub>O<sub>5</sub> at 20 K. Magnetic unit cell  $(a \times 2b \times c)$  is drawn.

Table 4. YBaFe<sub>2</sub>O<sub>5</sub> Structure Refinement Results from NPD Data

T (K):	20	280	300	320	340
λ (Å)	1.5402	2.0783	2.0783	2.0783	2.0783
$R_{ m wp}$	0.0671	0.0789	0.0907	0.0783	0.0748
$\chi^2$	4.78	1.38	1.91	1.25	1.17
space group	$Pmma^a$	$Pmma^a$	$Pmma^a$	$Pmmm^b$	$Pmmm^b$
a (Å)	8.0251(1)	8.0162(2)	8.0141(2)	3.93329(7)	3.93181(7)
b (Å)	3.83834(6)	3.85238(7)	3.85511(9)	3.91342(7)	3.91717(7)
c (Å)	7.5312(1)	7.5541(2)	7.5577(2)	7.5652(1)	7.5683(1)
$a_s-b_s(\text{Å})^c$	0.1742	0.1557	0.1520	0.0199	0.0146
$V(Å^3)$	231.983(8)	233.28(1)	233.49(1)	116.447(5)	116.563(5)
Fe(1) z	0.2542(4)	0.2568(9)	0.257(1)	0.2640(2)	0.2641(2)
Fe(2) z	0.2695(4)	0.2662(9)	0.265(1)		
O(1) z	0.003(1)	0.001(3)	0.000(3)	0	0
O(2a) z	0.3213(7)	0.324(2)	0.325(2)	0.3137(5)	0.3140(6)
O(2b) z	0.3132(7)	0.307(2)	0.308(2)		
O(3) x	0.0098(7)	0.012(1)	0.011(2)	0	0
O(3) z	0.3119(3)	0.3130(4)	0.3115(4)	0.3127(5)	0.3125(6)
$Y U_{iso} (\mathring{A}^2)$	0.0063(4)	0.0123(9)	0.013(1)	0.0152(8)	0.0137(7)
Ba $U_{iso}$ (Å <sup>2</sup> )	0.0039(6)	0.014(1)	0.017(1)	0.017(1)	0.0144(9)
Fe $U_{iso}/U_{eqv}$ $(A^2)^{d,e}$	0.0039(3)	0.0124(6)	0.0139(8)	0.0132	0.0122
$O(1)U_{\text{egv}} (A^2)^d$	0.0080	0.0138	0.0159	0.0181	0.0136
$O(2) U_{\text{egv}} $ $(\mathring{A}^2)^{af}$	0.0065	0.0123	0.0129	0.0186	0.0188
$O(3) U_{\text{egv}} $ $(\mathring{A}^2)^d$	0.0078	0.0128	0.0141	0.0146	0.0138
Fe $M_{\nu}$ ( $\mu_{\rm B}$ )	3.82(2)	3.41(3)	3.26(3)	2.88(2)	2.76(2)
Fe $M_z(\mu_B)$	0	0	0	-0.17(8)	-0.20(8)
Fe $M_{\text{Total}} (\mu_{\text{B}})$	3.82(2)	3.41(3)	3.26(3)	2.89(2)	2.77(2)

<sup>a</sup> Wyckoff positions for space group *Pmma* (nuclear cell) are: Ba at 2a (0,0,0); Y at 2c (0,0,\frac{1}{2}); Fe(1) and O(1) at 2f (\frac{1}{4},\frac{1}{2},\mathcal{Z}); Fe(2) at 2f (\frac{3}{4},\frac{1}{2},\mathcal{Z}); O(2a) at 2e (\frac{3}{4},0,\mathcal{Z}); O(2b) at 2e (\frac{1}{4},0,\mathcal{Z}); O(3) at 4j (x,\frac{1}{2},\mathcal{Z}). <sup>b</sup> Wyckoff positions for space group *Pmmm* (nuclear cell) are: Ba at 1a (0,0,0); Y at 1c (0,0,\frac{1}{2}); Fe at 2t (\frac{1}{2},\frac{1}{2},\mathcal{Z}); O(1) at 1f (\frac{1}{2},\frac{1}{2},0); O(2) at 2s (\frac{1}{2},0,\mathcal{Z}). <sup>c</sup> Orthorhombic distortion; refers to the single-perovskite-type subcell. <sup>d</sup>  $U_{eqv}$  values are given for those atoms where anisotropic displacement parameters were used in the refinement. The  $U_{eqv}$  values are defined as one-third the trace of the diagonal matrix describing the shape of the thermal ellipsoid. A complete list of the anisotropic displacement parameters is given in the Supporting Information. <sup>e</sup> The displacement parameters for Fe(1) and Fe(2) were constrained to be equal. An isotropic displacement parameter was used in the charge-ordered state, anisotropic displacement parameters were used for the MV state. <sup>f</sup> The displacement parameters for O(2a) and O(2b) were constrained to be equal.

## Caution !!!



### Things to keep in mind about powder diffraction measurements

Easy to do, but make sure you have a good powder!!!




FIG. 9. CCD image of Fe<sub>3</sub>Al powder (a) taken with stationary sample and (b) taken with sample rotation.


- Powder diffraction is excellent for getting the "big picture" but since intensities are spread over a sphere, small (but perhaps important) details are missed.
- Lose information about anisotropy

Physics 590 18


# **Thermodiffractometry**




## Detailed parameter dependence



High-flux mode leads to high data rates




**Figure 2.** DSC peaks for NdBaFe<sub>2</sub>O<sub>5</sub> (blue) and HoBaFe<sub>2</sub>O<sub>5</sub> (red) upon heating. Areas give latent heat  $(\Delta H)$  of the charge-ice melting.








# Single-crystal: more detail than powders Wide angle diffraction: Get an overview of everything















# Single crystal diffraction




## Triple-axis diffraction: focus in on specific points of interest



HB-1A 3-axis spectrometer





Orbital ordering in YVO<sub>3</sub>

## **Further references**



## General neutron scattering

- G. Squires, "Intro to theory of thermal neutron scattering", Dover, 1978.
- S. Lovesey, "Theory of neutron scattering from condensed matter", Oxford, 1984.
- R. Pynn, <a href="http://www.mrl.ucsb.edu/~pynn/">http://www.mrl.ucsb.edu/~pynn/</a>.

### Structural refinements

- GSAS <a href="http://www.ncnr.nist.gov/xtal/software/gsas.html">http://www.ncnr.nist.gov/xtal/software/gsas.html</a>
- FullProf http://www.ill.eu/sites/fullprof/

## How to get beam time

- Talk to one of us at Ames about your experiment
- We can identify a suitable instrument
- Talk to instrument scientist
- Write a beamtime request

Physics 590 23